

ON FARM FLEECE TESTING

On-farm **fleece testing** (PDF 1.41 MB) has become increasingly popular since 1999. On-farm testing may give some growers more flexibility in the way they manage their flocks. Fact sheets on the profitability of fleece testing, published by the Sheep CRC and Agriculture NSW in Australia, can be found [here](#). Since 2000, Australian Wool Innovation (AWI) and the Sheep CRC in Australia have invested heavily in R & D which has shown the ongoing financial benefits of individual animal measurements. Links to their sites are shown on our [Links](#) page.

[Info-bulletin 5.9](#) gives guidelines on the relative precision of different test methods. Refer to page 4 for information on fleece testing methods.

There is a range of options for on-farm fleece testing. Some growers have demonstrated that using this type of service close to shearing, primarily as a tool to assist with genetic selection, gives an added advantage that data has been gathered which allows them to split lines of wool based on micron and style grades. This can be quite simply done in the race in the yards a few hours (saves a muster) or days before shearing. A side benefit to this system is that the sheep can be shorn in micron groups relieving the pressure in the shed. Greasy fleece weights gathered at shearing can usefully add to the database.

SGS Wool Testing Services provides on-farm services using OFDA2000 instruments throughout New Zealand. The OFDA2000 instrument was developed as a fully-portable instrument, housed in a suitcase, specifically to allow diameter measurements to be made in real-time during animal selection or fleece classing see [Info-bulletin 3.4](#). In the simplest method of use, a single staple is selected from a consistent location on each fleece or animal, and this staple is prepared and measured in its greasy form. [Info-bulletin 3.7](#) addresses some of the issues relating to the measurement of greasy wool. The cycle of preparation and measurement takes 30 seconds or less, and throughputs of up to 1200 animals per day are commonly achieved in New Zealand by a single operator working alongside a selection team or shearing gang.

New Zealand growers have used the data from this service for both genetic selection and the objective preparation of lots for sale. Sale lots prepared in this manner are probably more uniform than lots prepared traditionally, and it remains to be seen whether buyers will pay a premium for this. Some growers have improved their wool return simply by breaking their clip into diameter-

based lines. [Info-bulletin 3.8](#) gives more detailed information on the relationship between fleece and core test diameter results. This technique will however only work when the market conditions are suitable and should not be relied upon as the principal benefit from the use of on-farm testing.

This service can also yield valuable data concerning the length-diameter profile of the wool, which is a very good indicator of the feed history of the sheep. [Click here](#) for further details on this and how such information can be used.

FOR ENQUIRIES

Email us at NZ.wool@sgs.com

Or contact us at:

48 Kemp Street, Kilbirnie

PO Box 15062

Wellington, New Zealand

Tel: +64.4.387.8565

Fax: +64.4.387.8651

© 2011 SGS. All rights reserved. The information contained herein is provided "as is" and SGS does not warrant that it will be error-free or will meet any particular criteria of performance or quality. Do not quote or refer any information herein without SGS' prior written consent. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.